

HOTLINE: Effective August 19, 2019

0050119

Lupus Comprehensive Reflexive Panel

LUPUS COMP

Note: Initial testing includes RF, C3, C4, and ANA. Specimens are screened for ANA using ELISA. If antibodies are detected, then an IFA titer will be added. If confirmed by IFA, then specimen will be tested for Thyroid Peroxidase (TPO) Antibody; Anti-Scl-70 (ENA), EIA; **Smith/RNP (ENA)** Antibody, IgG; **Smith (ENA)** Antibody, IgG; SSA 52 and 60 (Ro) (ENA) Antibodies, IgG; SSB (La) (ENA) Antibody, IgG; and Double-Stranded DNA (dsDNA) Antibody, IgG by ELISA; if Double-Stranded DNA (dsDNA) Antibody, IgG by ELISA result is Detected, then Double-Stranded DNA (dsDNA) Antibody, IgG by IFA (using *Crithidia luciliae*) is added. Additional charges apply.