
0050179 *Coccidioides* Antibody, IgG by ELISA**COCCI G****Specimen Required:** Collect: Serum Separator Tube (SST).Specimen Preparation: Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.3 mL)Storage/Transport Temperature: Refrigerated.Unacceptable Conditions: Contaminated, hemolyzed, or severely lipemic specimens.Stability (collection to initiation of testing): After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year
(avoid repeated freeze/thaw cycles)