
0050183 *Coccidioides immitis* Antibodies by Immunodiffusion**COCCL-PPT****Specimen Required:** Collect: Serum Separator Tube (SST).Specimen Preparation: Transfer 0.5 mL serum to an ARUP Standard Transport Tube (Min: 0.3 mL)Storage/Transport Temperature: Refrigerated.Unacceptable Conditions: Contaminated, hemolyzed, or severely lipemic specimens.Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year (avoid repeated freeze/thaw cycles)