

HOTLINE: Effective August 15, 2022

0050238 West Nile Virus Antibody, IgG by ELISA, CSF

WNIL IGG

Specimen Required: Collect: CSF.

Specimen Preparation: Transfer 2 mL CSF to an ARUP Standard Transport Tube. (Min: 0.3 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Bacterially contaminated, heat-inactivated, hemolyzed, or xanthochromic specimens.

Stability (collection to initiation of testing): Ambient: 8 hours; Refrigerated: 2 weeks; Frozen: 1 year (avoid repeated freeze/thaw cycles)