
0050341 **Immunoglobulin A, CSF****IGA CSF****Methodology:** Quantitative **Immunoturbidimetry****Specimen Required:** Collect: CSF.Specimen Preparation: Centrifuge and separate to remove cellular material. Transfer 1 mL CSF to an ARUP Standard Transport Tube. (Min: **0.5** mL)Storage/Transport Temperature: Refrigerated.Unacceptable Conditions: Grossly bloody or hemolyzed specimens.Stability (collection to initiation of testing): Ambient: **Unacceptable**; Refrigerated: **1 month**; Frozen: **6 months**