

HOTLINE: Effective February 16, 2021

0050356 IGM CSF Immunoglobulin M, CSF

Methodology: Quantitative Immunoturbidimetry

Performed: Wed, Sat Reported: 1-5 days

Specimen Required: Collect: CSF.

Specimen Preparation: Centrifuge and separate to remove cellular material. Transfer 1 mL CSF to an ARUP Standard Transport Tube.

(Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Grossly bloody or hemolyzed specimens.

Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 14 days; Frozen: 6 months