

HOTLINE: Effective September 3, 2019

0050375 Measles (Rubeola) Antibodies, IgG and IgM

MEASLE PAN

Reference Interval:

Test Number	Components	Reference Interval
0050380	Measles (Rubeola) Antibody, IgG	<p>Effective September 3, 2019</p> <p>13.4 AU/mL or less: Negative - No significant level of detectable measles (rubeola) IgG antibody. 13.5-16.4 AU/mL: Equivocal - Repeat testing in 10-14 days may be helpful. 16.5 AU/mL or greater: Positive - IgG antibody to measles (rubeola) detected, which may indicate a current or past exposure/immunization to measles (rubeola).</p>
0099597	Measles (Rubeola) Antibody, IgM	<p>0.79 AU or less: Negative - No significant level of IgM antibodies to measles (rubeola) virus detected. 0.80-1.20 AU: Equivocal - Repeat testing in 10-14 days may be helpful. 1.21 AU or greater: Positive - IgM antibodies to measles (rubeola) virus detected. Suggestive of current or recent infection or immunization. However, low levels of IgM antibodies may occasionally persist for more than 12 months post-infection or immunization.</p>