

HOTLINE: Effective September 3, 2019

**0050380 Measles (Rubeola) Antibody, IgG****MEASLES G****Reference Interval:**

Effective September 3, 2019

13.4 AU/mL or less	Negative - No significant level of detectable measles (rubeola) IgG antibody.
13.5-16.4 AU/mL	Equivocal - Repeat testing in 10-14 days may be helpful.
16.5 AU/mL or greater	Positive - IgG antibody to measles (rubeola) detected, which may indicate a current or past exposure/immunization to measles (rubeola).