
0050478**Rapid Plasma Reagin (RPR) with Reflex to Titer and TP-PA Confirmation****RPR PAN****Specimen Required:** Collect: Serum Separator TubeSpecimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.4 mL) Avoid freezing if possible.Storage/Transport Temperature: Refrigerated.Unacceptable Conditions: Plasma, CSF, or other body fluids.Stability (collection to initiation of testing): After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year (avoid repeated freeze/thaw cycles)