

0050577

Immunoglobulin G Subclasses (1, 2, 3, 4)

IGG SUB

Specimen Required: Collect: Serum separator tube.Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 2 mL serum to an ARUP Standard Transport Tube. (Min: 0.45 mL)Storage/Transport Temperature: Refrigerated.Stability (collection to initiation of testing): After separation from cells: Ambient: 2 hours; Refrigerated: 8 days; Frozen: 6 months