

0050627

Histoplasma Antibodies by CF & ID

HISTO PAN

Performed: Sun-Fri
Reported: 3-6 days

Specimen Required: Collect: Serum separator tube.

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL) Parallel testing is preferred and convalescent specimens **must** be received within 30 days from receipt of the acute specimens. **Mark specimens plainly as "acute" or "convalescent."**

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Contaminated or severely lipemic specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year (avoid repeated freeze/thaw cycles)