

---

**0050640 Protein Electrophoresis, Serum****SPEP****Performed:** Sun-Sat  
**Reported:** 1-3 days**Specimen Required:** Collect: Serum Separator Tube (SST).Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1.5 mL serum to an ARUP Standard Transport Tube. (Min: 1 mL)Storage/Transport Temperature: Refrigerated.Unacceptable Conditions: Plasma.Stability (collection to initiation of testing): After separation from cells: Ambient: **Unacceptable**; Refrigerated: 1 week; Frozen: 1 month