

HOTLINE: Effective August 19, 2019

0051495 Alpha Thalassemia (*HBA1* and *HBA2*) 7 Deletions

ALPHA THAL



Additional Technical Information



Out of Pocket Estimator

Out of Pocket Estimator

Performed: Varies
Reported: 7-10 days

Specimen Required: Collect: Lavender (EDTA), Pink (K₂EDTA), or Yellow (ACD Solution A or B).

Specimen Preparation: Transport 3 mL whole blood. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 6 months