

TEST CHANGE

Interleukin 12, Serum

0051530, IL12 DO

Specimen Requirements:	
Patient Preparation:	
Collect:	Serum separator tube, or plain red.
Specimen Preparation:	CRITICAL FROZEN: Separate serum from cells ASAP or within 2 hours of collection. Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.4 mL) <u>Freeze serum within 30 minutes.</u>
Transport Temperature:	CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered. Ship in an ARUP Standard Transport Tube.
Unacceptable Conditions:	Refrigerated specimens. Contaminated or heat-inactivated specimens.
Remarks:	
Stability:	After separation from cells: Ambient: 30 minutes; Refrigerated: Unacceptable; Frozen: 1 year
Methodology:	Quantitative Multiplex Bead Assay
Note:	Cytokine levels may demonstrate diurnal variation. For longitudinal comparison, it is recommended that cytokine levels be determined at the same time of day.

CPT Codes:	83520
------------	-------

New York DOH Approval Status:	This test is New York DOH approved.
-------------------------------	-------------------------------------

Interpretive Data:	▲
Results are used to understand the pathophysiology of immune, infectious, or inflammatory disorders, or may be used for research purposes.	

Reference Interval:	▲
Effective May 18, 2020 1.9 pg/mL or less Test Number	
Components	Reference Interval
Interleukin 12, Serum	2.6 pg/mL or less

Deleted Cells

Deleted Cells