

HOTLINE: Effective February 16, 2021

0051626 Epstein-Barr Virus Antibody to Viral Capsid Antigen, IgA EBV A

Reference Interval:

Effective February 16, 2021

Reference Interval	
0.7 U/L or less	Not Detected
0.8-1.2 U/L	Indeterminate – Repeat testing in 10-14 days may be helpful.
1.3 U/L or greater	Detected

HOTLINE NOTE: There is a unit of measure change associated with this test.

Change the unit of measure for component 0051626, EBV Antibody To Viral Capsid Antigen IgA from U/mL to U/L.