

HOTLINE: Effective September 7, 2021

0051668

Connective Tissue Diseases Profile

CONN

Methodology: Semi-Quantitative Enzyme-Linked Immunosorbent Assay
Semi-Quantitative Multiplex Bead Assay

Reference Interval: Effective May 18, 2015

Test Number	Components	Reference Interval		
0050085	Smith (ENA) Antibody, IgG			
		29 AU/mL or less	Negative	
		30-40 AU/mL	Equivocal	
		41 AU/mL or greater	Positive	
0050470	Smith/RNP (ENA) Antibody, IgG	Effective September 7, 2021		
		19 Units	Negative	
		20 to 39 Units	Weak Positive	
		40 to 80 Units	Moderate Positive	
		81 Units or greater	Strong Positive	
2012074	SSA 52 and 60 (Ro) (ENA) Antibodies, IgG			
		Test Number	Components	Reference Interval
			SSA-52 (Ro52) (ENA) Antibody, IgG	29 AU/mL or Less: Negative 30-40 AU/mL: Equivocal 41 AU/mL or greater: Positive
			SSA-60 (Ro60) (ENA) Antibody, IgG	29 AU/mL or Less: Negative 30-40 AU/mL: Equivocal 41 AU/mL or greater: Positive
0050692	SSB (La) (ENA) Antibody, IgG			
		29 AU/mL or less	Negative	
		30-40 AU/mL	Equivocal	
		41 AU/mL or greater	Positive	
0099592	Jo-1 Antibody, IgG			
		29 AU/mL or less	Negative	
		30-40 AU/mL	Equivocal	
		41 AU/mL or greater	Positive	
0099249	Ribosomal P Protein Antibody			
		29 AU/mL or less	Negative	
		30-40 AU/mL	Equivocal	
		41 AU/mL or greater	Positive	
0050714	Centromere Antibody, IgG			
		29 AU/mL or less	Negative	
		30-40 AU/mL	Equivocal	
		41 AU/mL or greater	Positive	
0050599	Scleroderma (Scl-70) (ENA) Antibody, IgG			
		29 AU/mL or less	Negative	
		30-40 AU/mL	Equivocal	
		41 AU/mL or greater	Positive	

HOTLINE NOTE: There is a unit of measure change associated with this test.
Change the unit of measure for component 0050470, Smith/RNP (ENA) Ab, IgG from AU/mL to Units.