
0051726*Leishmania* Antibody, IgG (Visceral Leishmaniasis)

LEISH IGG

Specimen Required: Collect: Serum separator tube.Specimen Preparation: Transfer 0.5 mL serum to an ARUP Standard Transport Tube. (Min: 0.1 mL)Storage/Transport Temperature: Refrigerated.Unacceptable Conditions: Serum containing glycerol or other viscous materials. Hemolyzed specimens.Stability (collection to initiation of testing): After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year