
0054441 Measles (Rubeola) Antibody, IgM, CSF**MEASLMCSF****Specimen Required:** Collect: CSF.Specimen Preparation: Transfer 0.5 mL CSF to an ARUP Standard Transport Tube. (Min: 0.2 mL)Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.Unacceptable Conditions: **Contaminated**, heat-inactivated, or hemolyzed specimens.Stability (collection to initiation of testing): Ambient: 8 hours; Refrigerated: 2 weeks; Frozen: 1 year