

HOTLINE: Effective April 3, 2019

0055284 Cysticercosis Antibody, IgG by ELISA

CYST SER

Performed: Sat
Reported: 1-8 days

Reference Interval:

Effective April 3, 2019

0.8 IV or less	Negative - No significant level of cysticercosis IgG antibody detected.
0.9 - 1.1 IV	Equivocal - Questionable presence of cysticercosis IgG antibody detected. Repeat testing in 10-14 days may be helpful.
1.2 IV or greater	Positive - IgG antibodies to cysticercosis detected, which may suggest current or past infection.

HOTLINE NOTE: There is a numeric map change associated with this test. There is also a unit of measure change associated with this test. Change the numeric map for component 0055284, Cysticercosis Ab, IgG by ELISA from XX.XX to XXX.X. Change the unit of measure for component 0055284, Cysticercosis Ab, IgG by ELISA from OD to IV.