

TEST CHANGE

Microsporidia Stain by Modified Trichrome

0060050, MICROST

Specimen Requirements:

Patient Preparation:

Collect: Stool. Recommended collection: 3 separate stool specimens within a 5-7-day period (an individual order must be submitted for each specimen).

Specimen Preparation: Preserve 2 g of stool within one hour of collection in AlcorFix (ARUP Supply #52059) available online through eSupply using ARUP Connect(TM) contact ARUP Client Services at (800) 522-2787. (Min: 1 g) Additional specimen collection instructions can be found at <https://www.aruplab.com/parasep>. Preserving in 10 percent formalin is also acceptable.

Transport Temperature: Room temperature.

Unacceptable Conditions: Unpreserved stool or specimens in any other preservative than indicated above.

Remarks:

Stability: Ambient: 9 months; Refrigerated: 9 months; Frozen: Unacceptable

Methodology: Qualitative Stain

Performed: ~~Mon-Fri~~ Sun-Sat

Reported: 1-~~3~~2 days

Note:

CPT Codes: 87207; 87015

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Refer to report.

Reference Interval:

Negative

