
0060101

Gram Stain

MS GRAM



Time Sensitive

Specimen Required: Collect: Any site or fluid.

Specimen Preparation: Transport labeled slide or specimen in Eswab transport media (ARUP Supply #45877) available online through eSupply using ARUP Connect™ or contact ARUP Client Services at (800) 522-2787.

Storage/Transport Temperature: Room temperature.

Remarks: Specimen source required.

Unacceptable Conditions: Blood, bone marrow, stool, vaginal specimens, or body fluid inoculated in culture bottles.

Stability (collection to initiation of testing): Ambient: 48 Hours; Refrigerated: 48 Hours; Frozen: Unacceptable