

HOTLINE: Effective August 16, 2021

0070103

Serum, C-Peptide

C PEP

Specimen Required: Patient Prep: Fasting specimen preferred.

Collect: Serum separator tube

Specimen Preparation: Allow specimen to clot completely at room temperature. Separate serum from cells ASAP. Submit specimen in an ARUP Standard Transport Tube.

Storage/Transport Temperature: Transport 1 mL serum, frozen. (Min: 0.5 mL)

Unacceptable Conditions: Grossly hemolyzed specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: 8 hours; Refrigerated: 7 days; Frozen: 1 month

Reference Interval:

Effective August 16, 2021

0.5-3.3 ng/mL

HOTLINE NOTE: There is a clinically significant charting name change associated with this test.

Change the charting name for component 0070103, C-Peptide, Serum or Plasma from C-Peptide, Serum or Plasma to Serum, C-Peptide.