
0080500

Alpha-1-Antitrypsin Phenotype (Includes Alpha-1-Antitrypsin)

A1A PHENO

Specimen Required: Collect: Serum Separator Tube (SST) or Plain Red.

Specimen Preparation: Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.3 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Grossly hemolyzed specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: 1 week; Refrigerated: 3 months; Frozen: 3 months (avoid repeated freeze/thaw cycles)