

---

**0080503****Lipoprotein Electrophoresis****LIPOELECT**

**Specimen Required:** Patient Prep: Patient should be fasting for 12-15 hours.

Collect: Serum separator tube or plain red.

Specimen Preparation: Transfer 1 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL) Store and ship refrigerated. DO NOT FREEZE.

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Heparin and frozen samples. Body fluids.

Stability (collection to initiation of testing): After separation from cells: Ambient: 24 hours; Refrigerated: 7 days; Frozen: Unacceptable