

0090299

Basement Membrane Zone and Cell Surface (Epithelial) Antibodies, IgG and IgA by IIF

EPITHELIAL



Immunodermatology Required Clinical
Information Form (Serum)

Methodology: Semi-Quantitative **Indirect Immunofluorescence (IIF)**

Specimen Required: Collect: Plain red or serum separator tube (SST).

Specimen Preparation: Transfer 2 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Hemolyzed or lipemic specimens. **Plasma.**

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 2 weeks; Frozen: Indefinitely

Interpretive Data: Refer to report