

TEST CHANGE

Heroin - Screen with Reflex to Confirmation/Quantitation - Serum or Plasma

0091203, HEROIN SP

Specimen Requirements:

Patient Preparation:

Collect: Plain ~~r~~Red or ~~gray~~ (sodium fluoride/potassium oxalate) ~~Gray~~
(Sodium Fluoride/Potassium Oxalate).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection.
Transfer 3 mL serum or plasma to an ARUP ~~standard transport~~
~~tube~~ ~~Standard Transport Tube~~. (Min: 1.0 mL)
Test is not performed at ARUP; separate specimens must be
submitted when multiple tests are ordered.

Transport Temperature: ~~Frozen~~ ~~Refrigerated~~. Also acceptable: ~~Refrigerated~~ ~~Frozen~~.

Unacceptable Conditions: ~~Polymer gel separation tube (SST or PST)~~ ~~Separator tubes~~.

Remarks:

Stability: Ambient: ~~2 days~~ ~~48 hours~~; Refrigerated: 2 weeks; Frozen: 1
month

Methodology: Qualitative Enzyme-Linked Immunosorbent Assay (~~ELISA~~) /
Quantitative High Performance Liquid Chromatography-
Tandem Mass Spectrometry

Note: If screen is positive, then confirmation will be added. Additional
charges apply.

CPT Codes: 80307; if reflexed, add 80356; 80361 (~~Reflexed Alt Code:~~
~~G0480~~)

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

~~Refer to~~ ~~By~~ report