
0091203

Heroin - Screen with Reflex to Confirmation/Quantitation - Serum or Plasma

HEROIN SP

Performed: Varies

Reported: 4-10 days

Specimen Required: Collect: Plain Red, Gray (Sodium Fluoride/Potassium Oxalate), Lavender (K₂ or K₃EDTA), or Pink (K₂EDTA).

Specimen Preparation: Transfer 3 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 1.3 mL)

Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.

Unacceptable Conditions: Separator tubes.

Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month