

HOTLINE: Effective August 19, 2019

0091224

LSD, Urine - Screen with Reflex to Confirmation/Quantitation

LSD URN

Performed: Varies

Reported: 8-18 days

Specimen Required: Collect: Random urine.

Specimen Preparation: Protect from light. Transfer 2 mL urine to an ARUP Amber Transport Tube. (Min: 0.9 mL)

Storage/Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Unacceptable Conditions: Specimens not protected from light. Specimens submitted in glass containers.

Stability (collection to initiation of testing): Ambient: 1 month; Refrigerated: 1 month; Frozen: 24 months