
0091224

LSD, Urine - Screen with Reflex to Confirmation/Quantitation

LSD URN

Performed: Varies
Reported: 8-18 days

Specimen Required: Collect: **Urine.**

Specimen Preparation: **Transfer** 2 mL urine to an **ARUP Standard Transport** Tube. (Min: **0.85 mL**)

Storage/Transport Temperature: Refrigerated. Also acceptable: **Frozen.**

Unacceptable Conditions: **Specimens** submitted in glass containers.

Stability (collection to initiation of testing): Ambient: **24 hours**; Refrigerated: 1 month; Frozen: **1 month**