
0091224

LSD, Urine - Screen with Reflex to Confirmation/Quantitation

LSD URN

Specimen Required: Collect: Urine.

Specimen Preparation: Transfer 2 mL urine to an ARUP Standard Transport Tube. (Min: 0.9 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.

Stability (collection to initiation of testing): Ambient: 24 hours; Refrigerated: 1 month; Frozen: 1 month