

TEST CHANGE

Fluvoxamine Quantitative, Serum or Plasma

0091234, FLUVOXAM

Specimen Requirements:

Patient Preparation:

Collect: Plain Red, Lavender (K2EDTA), or Pink (K2EDTA).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 3 mL serum or plasma to an ARUP Standard Transport Tube. (Min: 1.2 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Unacceptable Conditions: Separator tubes.

Remarks:

Stability: Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 6 months

Methodology: Quantitative Gas Chromatography

Performed: Varies

Reported: 5-8 days

Note:

CPT Codes: 80332 (Alt code: G0480)

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report