

HOTLINE: Effective January 1, 1900

0091585

Tin Total Quantitative, Serum or Plasma

TIN SP

Specimen Required: Collect: Royal Blue (K₂ EDTA), Royal Blue (Na₂ EDTA), or Royal Blue (No Additive).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an Acid Washed Transfer Vial (ARUP supply #54350) available online through eSupply using ARUP Connect™ or contact ARUP Client Services at (800) 522-2787. (Min: 0.4 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Unacceptable Conditions: Separator tubes.

Stability (collection to initiation of testing): Ambient: 28 days; Refrigerated: 28 days; Frozen: 28 days

CPT Code(s): 83789