

TEST CHANGE

Heroin - Screen with Reflex to Confirmation/Quantitation - Urine

0091586, HEROIN URN

Specimen Requirements:

Patient Preparation:

Collect: Random urine.

Specimen Preparation: Transfer 2 mL urine to an ARUP [standard transport tube](#)~~Standard Transport Tube~~. (Min: 1 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: [CRITICAL FROZEN](#)~~Frozen~~

Unacceptable Conditions: Thawed specimens.

Remarks:

Stability: Ambient: 24 hours; Refrigerated: [2 days](#)~~48 hours~~; Frozen: 3 months

Methodology: ~~Qualitative~~ Enzyme Immunoassay (EIA) / ~~Quantitative~~ High Performance Liquid Chromatography-Tandem Mass Spectrometry

Note: If screen is positive, then confirmation will be added. Additional charges apply.
~~Compare to Pain Management, Heroin Metabolite, Quantitative, with medMATCH, Urine and Pain Management, Heroin Metabolite, with Confirmation, with medMATCH, Urine~~

CPT Codes: 80307; if reflexed, add 80356; 80361 (~~Reflexed Alt Code: G0480~~)

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

[Refer to](#) ~~By~~ report