

---

**0091586**

**Heroin - Screen with Reflex to Confirmation/Quantitation - Urine**

**HEROIN URN**

**Specimen Required:** Collect: Random urine.

Specimen Preparation: Transfer 2 mL urine to an ARUP Standard Transport Tube. (Min: 1 mL)

**Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.**

Storage/Transport Temperature: Frozen

Unacceptable Conditions: Thawed specimens.

Stability (collection to initiation of testing): Ambient: **24 hours**; Refrigerated: 48 hours; Frozen: **3 months**