

TEST CHANGE

Pemphigoid Antibody Panel

0092001, PGOID PAN

Specimen Requirements:

Patient Preparation:

Collect: Plain red or serum separator tube.

Specimen Preparation: Transfer 2 mL serum to an ARUP [standard transport tube](#). ~~Standard Transport Tube~~. (Min: 0.5 mL)

Transport Temperature: Refrigerated/~~Ambient~~.

Unacceptable Conditions: Hemolyzed or lipemic specimens. Plasma.

Remarks: As a general rule, serum specimens should be shipped to the laboratory as soon as possible. Store refrigerated unless shipping promptly (within 2 hours). Transport at ambient temperature to arrive within 7 days. If 7-14 days until received in laboratory, store and ship refrigerated. If greater than 14 days, serum specimens must be stored and shipped frozen.

Stability: Ambient: 1 week; Refrigerated: 2 weeks; Frozen: Indefinitely

Methodology: Semi-Quantitative Indirect Immunofluorescence (IIF)/Semi-Quantitative Enzyme-Linked Immunosorbent Assay (ELISA)

Performed: Varies

Reported: ~~3~~4-9 days

Note: For specimens less than 0.5 mL, call the Immunodermatology Laboratory at [801-581-7139](#) ~~(866) 266-5699~~.

CPT Codes: 88346; 88350 x3; 83516 x2

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Refer to report.

Reference Interval:

By report

