

0092001

**Pemphigoid Antibody Panel**

**PGOID PAN**



Immunodermatology Required Clinical  
Information Form (Serum)

**Methodology:** Semi-Quantitative **Indirect Immunofluorescence (IIF)**/Semi-Quantitative Enzyme-Linked Immunosorbent Assay (**ELISA**)

**Specimen Required:** Collect: Plain red or serum separator tube.

Specimen Preparation: Transfer 2 mL serum to an ARUP Standard Transport Tube. (Min: 0.5 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Hemolyzed or lipemic specimens. **Plasma.**

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 2 weeks; Frozen: Indefinitely

**Interpretive Data:** Refer to report