

0092068

Hairstat 5 Reflexive Panel

HAIRSTAT 5

Specimen Required: Patient Prep: Ensure hair is not chemically treated or synthetic. Hair from the beard, underarms, chest, arms, legs or pubic hair may be collected. Body hair from different sites may be combined to get a final volume. Body hair and scalp hair should not be combined.

Collect: 100 mg hair. A kit must be ordered prior to collecting specimen (ARUP supply #40477). Available online through eSupply using ARUP Connect™ or by contacting Client Services at (800) 522-2787.

Specimen Preparation: Transport 100 mg hair (a ponytail 1.5 inches long and the diameter of a #2 pencil). (Min: 100 mg) **Specimen must have a tamper evident seal affixed. Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.**

Storage/Transport Temperature: Room temperature.

Unacceptable Conditions: **Specimens without a tamper evident seal affixed cannot be tested.**

Stability (collection to initiation of testing): Ambient: 1 year; Refrigerated: Undefined; Frozen: Undefined