
0092185

Drug Panel 7A, Urine - Screen with Reflex to Confirmation/Quantitation

CDASU 7A

Specimen Required: Collect: Random urine.

Specimen Preparation: Transfer 8 mL urine with no additives or preservatives to ARUP Standard Transport Tubes. (Min: 4 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Specimens exposed to repeated freeze/thaw cycles.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: **1 month**