

HOTLINE: Effective May 17, 2021

0092187

Drug Panel 9A, Urine - Screen with Reflex to Confirmation/Quantitation

CDASU 9A

Specimen Required: Collect: Random urine.

Specimen Preparation: Transfer 8 mL urine with no additives or preservatives in ARUP Standard Transport Tubes. (Min: 4 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Specimens exposed to repeated freeze/thaw cycles.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 1 month

HOTLINE NOTE: There is a reflexive pattern change associated with this test.

Add reflex to 3003726, Propoxyphene and Metabolite, Urine

Remove reflex from 2010468, Propoxyphene and Metabolite, Urine, Quantitative