

HOTLINE: Effective May 17, 2021

0092280

**Drugs of Abuse Test, Alcohol, Urine - Screen with Reflex to
Confirmation/Quantitation**

CDASUALC

Specimen Required: Collect: Random urine.

Specimen Preparation: Transport 10 mL urine. (Min: 2 mL)

Newborn minimum: 1 mL for the initial test and 1 mL for the confirmation of screen positive.

Storage/Transport Temperature: Room temperature.

Remarks: If applicable, indicate NEWBORN prominently on the test request form.

Unacceptable Conditions: Specimens exposed to repeated freeze/thaw cycles.

Stability (collection to initiation of testing): Ambient: 1 week; Refrigerated: 1 month; Frozen: 1 month