

TEST CHANGE

von Willebrand Factor Multimers

0092281, VWF MULTI

Specimen Requirements:

Patient Preparation:

Collect: Lt. blue (sodium citrate). ~~Special~~[Refer to Specimen Collection and Handling Hemostasis/Thrombosis Specimens guide located at https://www.aruplab.com/Specimen-Handling/SpecialSpecimenCollection/Hemostasis-Thrombosis.pdf for hemostasis/thrombosis specimen handling guidelines.](https://www.aruplab.com/Specimen-Handling/SpecialSpecimenCollection/Hemostasis-Thrombosis.pdf)

Specimen Preparation: Transfer 1 mL platelet-poor plasma to an ARUP ~~standard~~[transport tube](#). ~~Standard Transport Tube~~: (Min: 0.5 mL)

Transport Temperature: CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.

Unacceptable Conditions: Serum. EDTA plasma, clotted or hemolyzed specimens.

Remarks:

Stability: Ambient: 4 hours; Refrigerated: Unacceptable; Frozen ~~at -70 Degrees C: 6 months; Frozen at -20 Degrees C: 3 months~~

Methodology: Qualitative Electrophoresis

Performed: Mon-Fri

Reported: 4-11 days

Note:

CPT Codes: 85247

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

~~This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.~~

Reference Interval:

By report

Deleted Cells



*A nonprofit enterprise of the University of Utah
and its Department of Pathology*

Effective Date: **October 20, 2025**

