

**TEST CHANGE**

**Bullous Pemphigoid (BP180 and BP230) Antibodies, IgG by ELISA**

0092566, BP180 230G

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** Plain red or serum separator tube (SST).

**Specimen Preparation:** Transfer 2 mL serum to an ARUP [standard transport tube](#). ~~Standard Transport Tube~~. (Min: 0.5 mL)

**Transport Temperature:** Refrigerated/~~Ambient~~.

**Unacceptable Conditions:** Hemolyzed or lipemic specimens. Plasma.

**Remarks:** As a general rule, serum specimens should be shipped to the laboratory as soon as possible. Store refrigerated unless shipping promptly (within 2 hours). Transport at ambient temperature to arrive within 7 days. If 7-14 days until received in laboratory, store and ship refrigerated. If greater than 14 days, serum specimens must be stored and shipped frozen.

**Stability:** Ambient: 1 week; Refrigerated: 2 weeks; Frozen: Indefinitely

**Methodology:** Semi-Quantitative Enzyme-Linked Immunosorbent Assay ([ELISA](#))

**Performed:** Varies

**Reported:** 3-~~9~~<sup>11</sup> days

**Note:** For specimens less than 0.5 mL, call the Immunodermatology Laboratory at [801-581-7139](#)(~~866~~)-266-5699.

**CPT Codes:** 83516 x2

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

Refer to report.

**Reference Interval:**

By report



*A nonprofit enterprise of the University of Utah  
and its Department of Pathology*

Effective Date: **January 21, 2025**