

---

**0095200**

***Candida albicans* Antibodies IgA, IgG, and IgM by ELISA**

**CANDIDA AB**

**Specimen Required:** Collect: Serum separator tube.

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 0.5 mL serum to an ARUP Standard Transport Tube (Min: 0.15 mL)

Storage/Transport Temperature: Refrigerated

Unacceptable Conditions: Bacterially contaminated, heat-inactivated, hemolyzed, icteric, lipemic, or turbid specimens.

Stability (collection to initiation of testing): After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 6 months