

**TEST CHANGE** 

Arylsulfatase A, 24-Hour Urine

0095227, ARYLSULF A

Specimen Requirements:

**Patient Preparation:** 

Collect: 24-hour urine. Do not use preservatives. Keep refrigerated

during collection.

Specimen Preparation: From a well-mixed 24-hour collection, transfer Transfer 6 mL

urine to ARUP standard transport tubes. (Min. 2.5 mL) Test is

Effective Date: December 2, 2024

not performed at ARUP; separate specimens must be

submitted when multiple tests are ordered.

Transport Temperature: Refrigerated.

Unacceptable Conditions: Random urine.

Remarks: Indicate total volume.

Stability: Ambient: Unacceptable; Refrigerated: 2 weeks; Frozen:

Unacceptable

Methodology: Quantitative Colorimetry/Enzymatic Assay

Performed: Varies

Reported: 9-19 days

Note:

CPT Codes: 84311

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

By report