

TEST CHANGE

Antidepressant Panel Quantitative, Urine

0096048, ANTIDEP U

Specimen Requirements:

Patient Preparation:

Collect: ~~Urine~~ ~~Random urine~~.

Specimen Preparation: Transfer 2 mL urine to an ARUP standard transport tube. (Min: 0.7 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: ~~Frozen~~ ~~Refrigerated~~. Also acceptable: Room temperature or ~~refrigerated~~ ~~frozen~~.

Unacceptable Conditions:

Remarks:

Stability: Ambient: 1 week; Refrigerated: 11 days; Frozen: 2 weeks

Methodology: Quantitative Gas Chromatography / **Quantitative** Gas Chromatography-Mass Spectrometry (GC-MS)

Note: Panel includes: Amitriptyline, amoxapine, clomipramine, desmethylclomipramine, desipramine, doxepin, desmethyldoxepin, desmethyltrimipramine, fluoxetine, norfluoxetine, imipramine, maprotiline, mirtazapine, nortriptyline, protriptyline, trazodone, and trimipramine.

Desmethylsertraline (sertraline metabolite) and norcyclobenzaprine (cyclobenzaprine metabolite) are known interferences for protriptyline.

CPT Codes: 80332; 80337; 80338 (~~Alt code: G0480~~)

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval: