

---

0096657

Neutrophil Oxidative Burst Assay (DHR)

DHR

**Specimen Required:** Patient Prep: Collect control specimen from a healthy individual unrelated to patient at approximately the same time as and under similar conditions to the patient.

Collect: Green (Sodium or Lithium Heparin) (patient) AND Green (Sodium or Lithium Heparin) (control). **Patient and control specimens must be collected within 48 hours of test performance.**

Specimen Preparation: Transport 3 mL whole blood (patient) AND 3 mL whole blood (control) in original collection tubes. (Min: 1 mL (patient) AND 1 mL (control)) **Do not refrigerate or freeze. LIVE NEUTROPHILS REQUIRED.**

Storage/Transport Temperature: **CRITICAL ROOM TEMPERATURE.**

Stability (collection to initiation of testing): Ambient: 48 hours; Refrigerated: Unacceptable; Frozen: Unacceptable

**New York State Clients:** Ambient 24 hours; Refrigerated: Unacceptable; Frozen: Unacceptable

**HOTLINE NOTE:** Remove information found in the Unacceptable Conditions field.