

TEST CHANGE

Substance P, EIA

0098275, SUB P

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST).

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 5 mL serum to ARUP ~~standard transport tube~~ [Standard Transport Tubes](#) and freeze immediately. (Min: 1 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: CRITICAL FROZEN.

Unacceptable Conditions:

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 1 month

Methodology: Quantitative Enzyme Immunoassay [\(EIA\)](#)

Performed: Varies

Reported: ~~6-13~~[3-22](#) days

Note:

CPT Codes: 83520

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Reference Interval:

By report