

TEST CHANGE

17-Hydroxyprogesterone, Random Urine

0098378, 17OHPROGU

Specimen Requirements:

Patient Preparation: ~~Patient should not be on~~Discontinue any corticosteroid, ACTH, estrogen, or gonadotropin medications, if possible, for at least 48 hours prior to collection ~~of specimen.~~

Collect: ~~24-hour urine. No special preservatives required.~~
Urine

Specimen Preparation: ~~Mix specimen well.~~ Refrigerate during collection. Transfer 10 mL ~~aliquot of~~ urine to ARUP standard transport tubes. (Min: 5 mL) ~~Submit total volume.~~ Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: ~~Frozen. On dry ice is preferred. Separate specimens must be submitted when multiple tests are ordered.~~
Frozen.

Unacceptable Conditions:

Remarks:

Stability: Ambient: 1 hour; Refrigerated: 4 days; Frozen: 6 months

Methodology: Quantitative Radioimmunoassay (RIA)

Performed: Varies

Reported: 7-13 days

Note:

CPT Codes: 83498

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Reference Interval:

By report



*A nonprofit enterprise of the University of Utah
and its Department of Pathology*

Effective Date: **February 3, 2025**