

TEST CHANGE

Neuronal Cell Antibodies, <u>Quantitative,</u> CSF	
0098726, NEURON CSF	
Specimen Requirements:	
Patient Preparation:	
Collect:	<u>Cerebrospinal fluid (</u> CSF)-
Specimen Preparation:	Transfer 2 mL CSF to an ARUP <u>standard transport</u> <u>tubeStandard Transport Tube</u> . (Min: 1 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.
Transport Temperature:	Frozen. Also acceptable: Refrigerated.
Unacceptable Conditions:	
Remarks:	
Stability:	Ambient: Unacceptable; Refrigerated: 2 weeks; Frozen: Indefinitely
Methodology:	Quantitative Enzyme-Linked Immunosorbent Assay
Performed:	Varies
Reported:	3-8 days
Note:	
CPT Codes:	83520
New York DOH Approval Status:	Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.
Interpretive Data:	
Reference Interval:	
By report	
· ·	

HOTLINE NOTE: There is a component change associated with this test. One or more components have been added or removed. Refer to the Hotline Test Mix for interface build information.