

TEST CHANGE

Allergen, Food, Zucchini IgE

0099134, ZUCCHINI

Specimen Requirements:

Patient Preparation:

Collect: Plain red or serum separator tube (SST).

Specimen Preparation: Separate serum from cells ASAP or within 2 hours of collection. Transfer 0.5 mL serum plus 0.1 mL for each additional allergen ordered to an ARUP standard transport tube. Standard Transport Tube. (Min: 0.34 mL plus 0.04 mL for each allergen ordered)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Frozen. Room temperature. Also acceptable: Room temperature Refrigerated or refrigerated frozen.

Unacceptable Conditions: Hemolyzed, icteric, or lipemic specimens.

Remarks:

Stability: Ambient: 1 month; Refrigerated: 1 month; Frozen: 1 year

Methodology: Quantitative Enzyme Immunoassay (EIA)

Note:

CPT Codes: 86003

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Refer to By report